


Fast

Results in just a few minutes

- Results available in only 5-8 minutes
- Simple finger-prick (finger-stick) specimens
- Test is easily administered at point of care
- Each kit contains everything required to administer the test
- Each kit are individually packed, contains two test cassettes (IgG and IgM), dropper, dilution buffer, lancet and product manual.






Accurate

Dual IgG and IgM testing kits

- Separate IgG and IgM tests prevent a common problem in singlekit tests, where contaminants in the first test line can cause potential errors in the second test
- Positive test results indicate whether infection is in an early, strong or post-infection stage
- Clearly indicates if test returns an invalid result






Reliable

Conforms to global certification standards

- Developed by specialist labs in Taiwan
- Certified for use in many countries (CE Certification), with US FDA approval currently in process
- Storage requirements of 2–30°C make it easy to keep testing kits in optimum condition and fully effective
- Samples of each batch produced are kept by the factory for postproduction tracing of any issues that might warrant investigation







Value

Give you more than others

- Produced by the highest quality standard lab in Taiwan.
- Shelf life of 18 month enable piece of mind on inventory management.
- Testing kits can also be used to monitor Covid-19 vaccine effectiveness by measuring the presence of antibodies developed by the body's immune system after vaccination.







Service and Full Traceability

Each testing kit is traceable and accounted for

- Each testing kits are packaged individually with unique serial numbers, which make produce traceability extremely easy.
- In the unlikely event that the product has any kind of problem, the testing kit will be replace free of charge.

